Technical Specifications

Chambers	3
Compression Mode	4
Time Range	1min-24hours or continuous
Pressure Range	40-150mmHg
Hold Time	0-10s(Step by 1s)
Inteval Time	10-60s(Step by 1s)
Interface	5" touch screen
Rechargeable Battery	Yes (Lithium battery : 11.1V 9000mAh)
	Maximum Running Time: No Less than 8 Hours
Adaper	Yes
Input Voltage	AC100-240V 50-60Hz
Noise Level	<50 dB
Dimension	250(L)x184(W)x150(H)mm
Weight	2.5kg

Certification



CE (MDD)



Health Canada



Korea GMP



Guangzhou Longest Medical Technology Co., Ltd. Address: 301 & 401 of Building 2 & Building 3, No.96, Chuanggiang Road, Ningxi Street,

Address: 301 & 401 of Building 2 & Building 3, No.96, Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, Guangdong Province, 511399, China

Tel: 86-020-66353999

Fax: 86-020-66353920

E-mail: export@longest.cn **Website:** www.longestmedical.com





Health Life Longest Care






① Introduction of deep vein thrombosis (DVT)

Reason









Deep vein thrombosis is one of the most common complications of bed rest, prevention is urgently needed


Features

- 4 Compression Modes
- 5" Touch Screen
- Reusable & Disposable Sleeves
- O Gradient Pressure Setting
- Intermittent Pneumatic Compression
- Pressure of Every Chamber Adjustable
- Rechargeable Battery
- Powered by AC or DC



compression modes with reusable & disposable sleeves meet all your need

⊕Compression Modes



Ser Disposable Sleeve









Ser For Reusable Sleeve







