



Address: 301 & 401 of Building 2 & Building 3, No.96, Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, Guangdong Province, 511399, China

 $\textbf{E-mail:} \ export@longest.cn \qquad \textbf{Website:} \ www.longestmedical.com$ 









**Electrical Stimulation Machine LGT-2320BE** 

The LGT-2320BE is an advanced EMS machine, leveraging the power of low-frequency electrical stimulation therapy to achieve remarkable results in body sculpting. It not only efficiently contours and tones the body but also actively promotes muscle relaxation, relieving tension and soreness. Moreover, it plays a significant role in enhancing skin elasticity, rejuvenating the skin's appearance and suppleness.





#### **Body Sculpting**

It electrically stimulates the muscles of the abdomen, buttocks, thighs, and other areas to promote fat metabolism and muscle toning, thereby creating a beautiful body curve



#### **Facial Rejuvenation**

It stimulates the contraction of facial muscles, lifts and tightens the skin, reduces wrinkles, and promotes blood circulation and collagen production, enabling the face to regain its glow.



#### **Muscle Building**

Whether for post-injury or illness rehabilitation or fitness requirements, it accurately stimulates the muscles to strengthen muscle strength and endurance.



#### **Postpartum Recovery**

It assists new mothers in restoring the function of muscles damaged during pregnancy and childbirth and simultaneously alleviates problems such as abdominal looseness, expediting the rapid recovery of body functions.

# **Benefits of Incorporating EMS Machine for Your Beauty Treatment Services**

- No Needles
- No Anesthesia
- Safe and Reliable
- Enhanced Client Comfort
- Improved Treatment Efficiency
- Cost-effective Investment
- Versatile Application
- Time-efficient
- Visible and Lasting Results

### **Features**



## **Accessories**



Five kinds of pen electrodes for detailed facial treatment.



Twelve pieces of silicon electrodes applicable to different treatment sites.



**Ball Point Electrode** 



\*Spherical Electrode




\*Plate Electrode



\*Y-shape Electrode



\*Dual-PointBall Electrode



\*Glove Electrode(1pair)



**Silicon Electrode** 



**Self-Adhesive Electrode** 

<sup>\*</sup> Optional Accessories